

## **Product Information**


### **0.2-10ul Micro Tips, Medi**

#### **Description**

Biobasic Inc. made pipette tips by high purity virgin polypropylene and are designed for sensitive applications. They are compatible with robotic and manual pipetting systems and are tested for compatibility with Eppendorf®, Gilson®, and Biobasic pipettes.

#### **Features**

- Ideal for Precise, Accurate Pipetting: Tips are straight and feature clean, Uniform bores and hydrophobic surfaces to reduce sample retention and increase accuracy and reproducibility
- Available in variety of sizes.
- Autoclavable.



Item No.	Description	Package
BT216	Natural/bag	1000/bag
BT216-S	Natural/bag/sterile	100/bagx10
BT216-RS	Natural/rack/sterile Certificated RNase, DNase-Free.	96/rack

#### **Technical Specifications**

Item No.	BT216	BT216-S	BT216-RS
Package	1000/Bag	100/Bagx10	96/RACK
Volume	0.2 -10µl	0.2-10µl	0.2-10µl
Graduations	NO	NO	NO
Color	Natural	Natural	Natural
Plastic materials	Polypropylene	Polypropylene	Polypropylene
Sterile	NO	YES	YES
DNase, RNase Free	NO	YES	YES
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C
Length	36 mm	36 mm	36 mm
Storage condition	Room temperature	Room temperature	Room temperature

#### **Comparison of Related Tips**



QF 24 Rev 1