

Product information

QF 24 V4
 V1 July 2021

Compact Balances

Catalog #: BCM1300/1301/1302/1303/1304/1305
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Compact balances for a wide range of general weighing applications. Six models cover up to 10000 grams. Multiple program options, including standard weighing, count mode, percent, etc.

Features:

- RS232 port for data output
- Gold-plated ceramic capacitance sensor for accurate weighing
- Enhanced shock resistance and overload protection
- Large, user-friendly interface

Technical Data:

	BCM1300	BCM1301	BCM1302	BCM1303	BCM1304	BCM1305
Capacity	120 grams	300 grams	500 grams	1200 grams	5000 grams	10000 grams
Readability	0.01 grams			0.1 grams		
Pan Size (round)	11.5 cm / 4.5"			16 cm / 6.3"	18 cm / 7"	
Draft Shield	Yes, plastic ring			No		
Weighing Units	g, oz, ct					
Calibration	External					
Display	Backlit LCD					
Dimensions	18.8 x 28 x 8 cm / 7.1 x 11 x 3.1 inches					

Ordering Information:

BCM1300 - Compact Balance, 120 grams, 115V

BCM1301 - Compact Balance, 300 grams, 115V

BCM1302 - Compact Balance, 500 grams, 115V

BCM1303 - Compact Balance, 1200 grams, 115V

BCM1304 - Compact Balance, 5000 grams, 115V

BCM1305 - Compact Balance, 10000 grams, 115V



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.