

## Product information

QF 24 V4  
 V1 July 2021

# CO<sub>2</sub> Orbital Shaker

**Catalog #:** BCM1625  
**Size:** 1 UNIT  
**Storage:** 18-25°C

### Product Description:

Orbital shaker with flat mat, suitable for use in CO<sub>2</sub> incubators and other extreme environments. Ideal for aeration and mixing. Includes remote control. Speed can be adjusted. 2 year warranty.

### Features:

- Suitable for use in CO<sub>2</sub> incubators
- Includes remote control with 60.5"/154 cm cord length
- Adjustable speed
- Includes non-slip rubber mat platform
- Safe for incubator or cold room use
- Magnetic clamp platform for flasks and tube racks is not included with BCM1625



**BCM1625** - Shaker CO<sub>2</sub> with remote controller and rubber mat platform, 13" x 12", 115V

### Technical Data:

<b>Speed Range</b>	30-300rpm; adjustable in increments of 1rpm
<b>Orbit</b>	19mm
<b>Load Capacity</b>	5 x 1L or 14lbs (6.4kgs)
<b>Shaker Dimensions (WxDxH)</b>	12 x 13.4 x 3.5 inches 30 x 34 x 9 centimeters
<b>Remote Dimensions (WxDxH)</b>	7.2 x 6.7 x 3.5 inches 18 x 17 x 9 centimeters
<b>Platform Size</b>	13 x 12 inches 33 x 30 centimeters
<b>Weight</b>	25 lbs / 11.5 kg
<b>Timer</b>	1min to 48hr/continuous
<b>Electrical</b>	120V, 50-60Hz, 600W
<b>Operating Temperature</b>	0°C to +65°C
<b>Warranty</b>	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
 NOT INTENDED FOR HUMAN OR ANIMAL USE.