

Product information

QF 24 V4
 V1 July 2021

Shaking Water Bath

Catalog #: BCM1404
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Shaking water bath with universal spring platform for tubes, bottles and flasks. Hinged, removable lid has a gabled design to direct condensation away from center. 2 year warranty.

Features:

- Universal spring platform included
- Seamless stainless steel chamber
- Digital temperature, speed and timer controls



BCM1404 - Shaking Water Bath

Technical Data:

Temperature Range	+5°C to +80°C, 0.1°C incr.
Temperature Accuracy	±0.3°C (at 37°C)
Temperature Uniformity	±0.2°C (at 37°C)
Speed Range	30 – 200 rpm (1rpm incr.)
Shaking Orbit	25mm
Platform	Universal Spring Platform
Exterior Dimensions (WxDxH)	12.2 x 13.25 x 12 inches 32 x 34 x 30 cm
Timer	1min to 48hr / Continuous
Weight	31 lbs / 14 kg
Electrical	115V, 50-60 Hz, 300W
Warranty	2 years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.