

Product information

QF 24 V4
 V1 July 2021

Digital Dry Bath

Catalog #: BCM1411
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Dry bath with advanced microprocessor controls, timed operation and a removable hinged lid. The bath contains two isolated, independently controlled chambers allowing the bath to serve as two baths in one unit. Settings for either side may be adjusted at any point in operation without affecting the settings of the other side. Dividing wall acts as a barrier to prevent air transfer between sides. 2 year warranty.

Features:

- Hinged cover
- Two isolated chambers with independent temperature and timer controls
- Digital controls



BCM1411 - Digital Dry Bath

Technical Data:

Temperature Range	Ambient +5°C to +105°C Increments of 0.1°C Accuracy of ±0.1°C
Temperature Uniformity	±0.2°C
Timer	1min to 99hr 59min
Dimensions (WxDxH)	6.3 x 14.5 x 5 inches 16 x 36.5 x 13 centimeters
Electrical	115V, 50-60Hz, 240W
Operating Temperature	+4°C to +65°C
Warranty	2 years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.