

## Product information

QF 24 V4  
 V1 July 2021

# Autoclave

**Catalog #:** BCM1541/1542  
**Size:** 1 UNIT  
**Storage:** 18-25°C

### Product Description:

Research autoclave, with 16L interior. Includes mechanical and electrical lock to prevent accidental opening during operation. BCM1542 is an optional mini printer that provides a record of sterilization parameters. 1 year warranty.

### Features:

- Fully automatic – 3 available programs
  - Sterilization of liquids, unwrapped instrument, and wrapped instrument settings
- All settings are adjustable
- Mechanical and electrical safety locks prevent door from opening until interior reaches zero PSI
- Built-in water tank

### Technical Data:

<b>Chamber Volume</b>	16L
<b>Max. Pressure</b>	29 PSI / 2 Bar
<b>Weight</b>	89lbs/40kg
<b>Dimensions (WxDxH)</b>	22 x 17.5 x 15.7 inches 56 x 44.5 x 40 centimeters
<b>Chamber Dimensions</b>	9 x 13.5 inches 23 x 35 centimeters
<b>Tray Dimensions</b>	6.5 x 10.5 inches 16.5 x 27 centimeters
<b>Electrical</b>	115V, 50-60Hz, 1600W
<b>Warranty</b>	1 Year



### Ordering Information:

**BCM1541** - Research Autoclave, 16L, 115V

**BCM1542** - Printer for BCM1541, 115V



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
 NOT INTENDED FOR HUMAN OR ANIMAL USE.