

Product information

QF 24 V4
 V1 July 2021

Adjustable Platform Rocker

Catalog #: BCM1614
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Platform rocker with flat mat. Speed and tilt can be adjusted. 2 year warranty.

Features:

- Adjustable tilt and speed
- Large work surface
- Safe for incubator or cold room use



BCM1614 - Platform Rocker, fixed speed,
 12" x 12" with flat mat, 115V

Technical Data:

Speed	2-30rpm
Tilt Angle	0° to ±30°
Load Capacity	8 lbs / 3.6 kg
Dimensions (WxDxH)	14 x 12 x 8 inches 35 x 30 x 20 centimeters
Weight	7.5 lbs / 3.6 kg
Electrical	115V, 50-60Hz, 50W
Operating Temp Range	+4°C to +65°C
Platform Size	14 x 12 inches / 35 x 30 cm
Warranty	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.