

Product information

QF 24 V4
 V1 July 2021

Variable 3D Nutating Shaker

Catalog #: BCM1617
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Multifunctional mixer that provides a 3D nutating motion for thorough mixing. Speed and angle are adjustable. Instantly adjust platform tilt angle by gripping platform with two hands and moving platform to the desired position. Brushless, maintenance-free motors. Includes both flat and dimpled mats. 2 year warranty.

Features:

- Includes dimpled mat for use with centrifuge or blood tubes in addition to flat mat
- Instantly adjust platform tilt angle
- Mat is autoclavable
- Compatible with BCM1645 bungee cords to secure trays to platform
- Safe for incubator or cold room use



BCM1617 - Mixer 3D Nutating Shaker, fixed speed, 10" x 7.5" platform, dimpled mat, 115V

Technical Data:

Speed	2-30rpm
Tilt Angle	0 to $\pm 30^\circ$
Load Capacity	5 lbs / 2.3 kg
Weight	7.5 lbs / 3.4 kg
Electrical	115V, 50-60Hz, 50W
Dimensions (WxDxH)	12 x 12 x 8 inches 30 x 30 x 20 centimeters
Platform Size	12 x 12 inches 30 x 30 centimeters
Operating Temperature	+4°C to +65°C
Warranty	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.