

## Product information

QF 24 V4  
 V1 Sept 2021

# 3D Blot Rocker

**Catalog #:** BCM1619  
**Size:** 1 UNIT  
**Storage:** 18-25°C

### Product Description:

The 3D Blot Rocker provides the perfect 3-dimensional motion for all blot related applications. Its tilt angle and speed have been optimally set for gentle, but thorough mixing gel trays, boxes and other flat vessels. The 3D Blot Rocker comes standard with autoclavable flat mats and is ideal for blot washes and gel/staining. The 3D Blot Rocker has been designed to conserve bench space. The 3D Blot Rocker accepts stackable platforms and handle loads up to 2kg. It has maintenance free brushless motors and is safe for use in cold rooms and incubators.

### Features:

- Flat mat is autoclavable
- Gentle three-dimensional rotations thoroughly rinse gels
- Compatible with BCM1645 bungee cords to secure trays to platform
- Safe for incubator or cold room use

### Technical Data:

<b>Speed</b>	12rpm
<b>Tilt Angle</b>	±8°
<b>Load Capacity</b>	3.5lbs / 1.6 kg
<b>Weight</b>	4.4 lbs / 2 kg
<b>Electrical</b>	115V, 50-60Hz, 50W
<b>Dimensions (WxDxH)</b>	12 x 12 x 8.3 inches 30 x 30 x 20.4 centimeters
<b>Platform Size</b>	12 x 12 inches 30 x 30 cm
<b>Operating Temperature</b>	+4°C to +50°C
<b>Warranty</b>	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
 NOT INTENDED FOR HUMAN OR ANIMAL USE.