



CERTIFICATE OF ANALYSIS

Product	PROTEASE INHIBITOR COCKTAIL (For Bacteria Cell, EDTA Free, 1ml)
Product Code	BS380
Lot #	Q53080R0X

Components:

Component	BS380
PMSF	100mM
Bestatin	2 mM
Pepstatin A	0.3 mM
E-64	0.3 mM

Storage:

Store cocktail at -20°C . Cocktail is dissolved in DMSO containing little deionized water. The cocktail may form crystals. If there are crystals, place it at RT for 5 minutes.

This is to certify that units of the above mentioned lot number were tested and found to comply with the specifications of the grade listed. Bio Basic expressly disclaims all warranties, expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. The following are the actual analytical results obtained:

Result name	Units	Specification	Test Value
Functional Test	Pass/Fail	=Pass Test	Passed