



CERTIFICATE OF ANALYSIS

Product	PROTEASE INHIBITOR COCKTAIL (For Bacteria Cell, with EDTA, 1ml)
Product Code	BS381
Lot #	L120ROJ

Components:

Component	BS381
PMSF	100 mM
EDTA	100 mM
Bestatin	2 mM
Pepstatin A	0.3 mM
E-64	0.3 mM

This is to certify that units of the above mentioned lot number were tested and found to comply with the specifications of the grade listed. Bio Basic expressly disclaims all warranties, expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. The following are the actual analytical results obtained:

Result name	Units	Specification	Test Value
Appearance		Report	Kit
Function Test	Pass/Fail	=Pass Test	Pass test

Storage:

Store cocktail at -20°C . Cocktail is dissolved in DMSO containing little deionized water. The cocktail may form crystals. If there are crystals, place it at RT for 5 minutes.