



## CERTIFICATE OF ANALYSIS

Product	PROTEASE INHIBITOR COCKTAIL (For Mammalian Cell, 1ml)
Product Code	BS386
Lot #	P80486R0T
Retest Date	February, 2023

### Components:

Component	BS386
PMSF	100 mM
Bestatin	1 mM
Pepstatin A	1.5 mM
E-64	1.4 mM
Aprotinin	0.08 mM
Leupeptin	1 mM


This is to certify that units of the above mentioned lot number were tested and found to comply with the specifications of the grade listed. Bio Basic expressly disclaims all warranties, expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. The following are the actual analytical results obtained:

Result name	Units	Specification	Test Value
Appearance		Report	Kit
Functional Test	Pass/Fail	=Pass Test	Pass test

### Storage:

Store cocktail at  $-20^{\circ}\text{C}$ . Cocktail is dissolved in DMSO containing little deionized water. The cocktail may form crystals. If there are crystals, place it at RT for 5 minutes.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:   
Title: General Manager  
Date: April 26, 2022

QF 24 Rev 1