



CERTIFICATE OF ANALYSIS

| | |
|--------------|---------------------------------------|
| Product | PBS Blocking buffer with non fat milk |
| Product Code | SB0625 |
| Lot No | N/A |

| Test Items | Specifications | Actual Results |
|---|----------------|----------------|
| Appearance | Powder | |
| Component of per PK | | |
| Sodium phosphate, dibasic, heptahydrate | 2.17G | |
| Potassium chloride | 0.2G | |
| Sodium chloride | 8G | |
| Potassium phosphate, monobasic, Anhydrous | 0.2G | |
| Non-Fat Powdered Milk | 50G | |

Storage: -15 to -20°C