



Certificate of Analysis

| | |
|--------------|------------------------------------|
| Product | Silica gel with moisture indicator |
| Product Code | DSG110 |
| CAS# | 112926-00-8 |
| Lot No | Q4024860 |
| Retest Date | January 04, 2026 |

| Test Items | Specifications | Results |
|------------------------------|---------------------------------------|-------------|
| Appearance | Blue to light blue Spherical particle | Conforms |
| Rate of water | ≤10.0% | 9.3% |
| Loss on drying | ≤5.0% | 3.1% |
| Spherical particle Pass rate | ≥82% | 82% |
| Adsorption | | |
| RH=20% | ≥8% | 8.7% |
| RH=35% | ≥13% | 13.7% |
| RH=50% | ≥ 20% | 21.4% |
| Colour | | |
| RH=20% | Blue or light blue | Blue |
| RH=35% | Purple or Amaranthine | Amaranthine |
| RH=50% | Light red | Light red |

Storage: 18~25°C. Protect from moisture

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: April 11, 2022