

## Technical bulletin

QF 24 V4  
V2 May 2023

# Lysis Buffers Usage

SKU	Description	Size	Use
BSP001.SIZE.50Preps	Extract-EZ G, Nyclear & Cytoplasmic Protein Extracton Kit	50 Preps	Mammalian Cells or Tissues
BSP002.SIZE.50Preps	Membrane, Nuclear & Cytoplasmic Protein Extraction Kit	50 Preps	Mammalian Cells or Tissues
BSP003.SIZE.50Preps	Total Protein Extraction Kit	50 Preps	Mammalian Cells or Tissues
BSP004.SIZE.20Preps	Plant Protein Extraction Kit	20 Preps	Plant Tissue
BSP005.SIZE.50Preps	Extract-EZ K1, membrane and Cytoplasmic Protein Extraction Kit	50 Preps	Mammalian Cells or Tissues
BSP006-20ml.SIZE.20ml	Tissue Proteins Lysis Buffer, 1X Strength Solution	20 ml	Mammalian Tissue
BSP006-5x20ml.SIZE.5x20ml	Tissue Proteins Lysis Buffer, 1X Strength Solution	5x20 ml	Mammalian Tissue
BSP009.SIZE.50Preps	Nuclear Protein Fractionation Kit	50 Preps	Mammalian Cells or Tissues
BSP011.SIZE.50Preps	TCA Precipitating Protein Kit	50 Preps	Protein Extraction from Animal Tissues, Cells or Protoplasts
BSP012.SIZE.100Preps	Universal Protein Precipitating Reagent	100 Preps	Protein Extraction from Animal Tissues, Cells or Protoplasts
BSP013.SIZE.1kit	Extract-EZ C Yeast Protein Extraction Kit	1 kit	Yeast Sample
BSP022-20ml.SIZE.20ml	Animal Cells Protein Lysis Buffer, 1x Strength Solution	20 ml	Animal Cell
BSP022-5x20ml.SIZE.5x20ml	Animal Cells Protein Lysis Buffer, 1x Strength Solution	5x20 ml	Animal Cell
BSP023-25ml.SIZE.25ml	Bacterial Protein Lysis Buffer, 1x Stength Solution	25 ml	Bacterial Cell
BSP023-4x25ml.SIZE.4x25ml	Bacterial Protein Lysis Buffer, 1x Stength Solution	4x25 ml	Bacterial Cell
BSP026-20ml.SIZE.20ml	Yeast Protein Lysis Buffer, 1x Strength Solution	20 ml	Yeast Sample
BSP026-5x20ml.SIZE.5x20ml	Yeast Protein Lysis Buffer, 1x Strength Solution	5x20 ml	Yeast Sample
BSP051.SIZE.50Preps	Extract-EZ L1, Cytosol & Mitochondria Protein Extraction Kit	50 Preps	Mammalian Cells or Tissues
BSP055.SIZE.100Assays	SDS Dectection Kit	100 Assays	Animal Samples & other Samples involving SDS Contamination
BSP058.SIZE.50Preps	FFPE Protein Extraction Kit	50 Preps	Paraffin Embedded Tissue
BSP061.SIZE.2000Assays	Phosphoprotein Phosphate Estimation Assay Kit	2000 Assays	Protein Extraction from Animal Tissues, Cells or Protoplasts
BSP073.SIZE.50Preps	Extract-EZ E1, Total Subcellular Protein Extraction Kit	50 Preps	Animal Cell
PL005-10ml.SIZE.10ml	RIPA Lysis & Extraction Buffer I, 1x Strength Solution	10 ml	Mammalian Cells or Tissues
PL005-5X10ml.SIZE.5x10ml	RIPA Lysis & Extraction Buffer I, 1x Strength Solution	5x10 ml	Mammalian Cells or Tissues
PL006-10ml.SIZE.10ml	RIPA Lysis & Extraction Buffer II, 1x Strength Solution	10 ml	Mammalian Cells or Tissues
PL006-5x10ml.SIZE.5x10ml	RIPA Lysis & Extraction Buffer II, 1x Strength Solution	5x10 ml	Mammalian Cells or Tissues
PL007-10ml.SIZE.10ml	RIPA Lysis & Extraction Buffer III, 1x Strength Solution	10 ml	Mammalian Cells or Tissues
PL007-5x10ml.SIZE.5x10ml	RIPA Lysis & Extraction Buffer III, 1x Strength Solution	5x10 ml	Mammalian Cells or Tissues
PL008-10ml.SIZE.10ml	RIPA Lysis & Extraction Buffer IV, 1x Strength Solution	10 ml	Mammalian Cells or Tissues
PL008-5x10ml.SIZE.5x10ml	RIPA Lysis & Extraction Buffer IV, 1x Strength Solution	5x10 ml	Mammalian Cells or Tissues