



Four Source Power Supply

(Cat. #EDK621)



Power any combination of Electrophoresis Apparatus units with this mighty 750mA power supply! Features an easy-to-use, fully programmable interface for setting voltage, current or timer control with each parameter displayed in real-time. Programs may be paused or resumed at any point. Run experiments in the least time possible with this powerful and versatile unit!



10-300 V for **1 to 4 units**



Cat. #EDK621

Features:

- Max Voltage: 300 Volts
- Voltage Steps: 10-300 Volts
- Max Current: 750 Milliamps
- Output Type: Variable Voltage
- Lead Inputs: 4 Sets, Recessed, Color Coded
- Fuse: 3.15 Amp Glass Fuse
- Input Power: 60 Hz, 110 V

