



## Product Information Mini Digital Dry Bath

### Product information for BCM1405/BCM1406:

#### Introduction

Personal benchtop digital dry bath. Dividing wall acts as a barrier to prevent air transfer between sides. Simple, accurate “set-and-forget” controls. Bath and lid temperature for BCM1406 may be selected independently. No timer controls. 2 year warranty.

**Please note, BCM1405 and BCM1406 do not come with blocks included. See below for product codes of blocks that may be purchased.**



\* Heat block sold separately.

#### Features:

- BCM1406 features a heated lid to prevent condensation and evaporation
- Less than five inches wide, fits on the most crowded benchtop
- Digital controls

#### Ordering Information

| Product Code | Description                                       |
|--------------|---|
| BCM1405      | Mini dry bath, 100-240V (US Plug)                 |
| BCM1406      | Mini dry bath with heated lid, 100-240V (US Plug) |

#### Technical Data:

|                               |   |                           |  |
|-------------------------------|---|---------------------------|--|
| <b>Temperature Range</b>      | Ambient +5°C to +100°C<br>Increments of 0.1°C<br>Accuracy of ±0.5°C | <b>Dimensions (WxDxH)</b> | 4 x 6 x 4.5 inches<br>10.5 x 55 x 11.5 centimeters |
| <b>Temperature Uniformity</b> | ±0.2°C  | <b>Weight</b>             | 2lbs / 1kg   |
| <b>Electrical</b>             | Universal voltage 100-240V, 50-60Hz, 50W                            | <b>Warranty</b>           | 2 Years  |



## Ordering Information for Blocks:

| Product Code   | Description  |
|----------------|--|
| BCM1405-BCU    | Block. 1.25in/32mm in height. Custom drilled with up to 20 holes |
| BCM1405-BCV    | Block, for 8 cuvettes (12.5x12.5x32mm)                           |
| BCM1405-B02    | Block, 40 x 0.2ml tubes or 5 PCR strips of 8 tubes each          |
| BCM1405-B05    | Block, 24 x 0.5ml centrifuge tubes                               |
| BCM1405-B1213  | Block, 15 x 12-13mm tubes, 60mm tall, no lid                     |
| BCM1405-B1232  | Block, 15 x HPLC or Cryovial tubes                               |
| BCM1405-B15    | Block, 15 x 1.5ml centrifuge tubes                               |
| BCM1405-B1520  | Block, 15 x 1.5ml or 2.0ml centrifuge tubes (conical)            |
| BCM1405-B15HL* | Block, 15 x 1.5ml Centrifuge Tubes for Dry Bath with Heated Lid* |
| BCM1405-B5MT   | Block, 6 x 5ml centrifuge tubes (17mm diameter)                  |
| BCM1405-B150   | Block, 4 x 15ml centrifuge tubes                                 |
| BCM1405-B500   | Block, 2 x 50ml centrifuge tubes                                 |

\* BCM1405-B15HL is for use with BCM1406 only.