



Product Information

Mini Incubator

Product information for BCM1415:

Introduction

Heats or cools samples. Digital temperature control allows one to “set-and-forget”, eliminating the need for fine-tuning. Includes two adjustable shelves. 2 year warranty.

Features:

- Digital temperature control
- Compact, personal-sized footprint
- Accepts bottles and flasks up to 2L



Ordering Information

Product Code	Description
BCM1415	Mini digital incubator with heating and cooling, 115V

Technical Data:

Temperature Uniformity	$\pm 1^{\circ}\text{C}$ (at 37°C)	Interior Dimensions (WxDxH)	10.3 x 9.3 x 12.8 inches 26 x 23.5 x 33 centimeters
Temperature Range	Ambient -15°C to 60°C , 1°C increments	Exterior Dimensions (WxDxH)	13.2 x 14.5 x 18.7 inches 33.5 x 37 x 48 centimeters
Weight	15lbs/6.5kg	Warranty	2 Years
Electrical	115V, 50-60Hz		