



## **Product Information**

### **Mini Blot 3D Rocker**

#### ***Product information for BCM1618:***

#### **Introduction**

Compact fixed-speed 3D rocker for blot washing. Optional accessories are available. 2 year warranty.

#### **Features:**

- Compact
- Gentle three-dimensional rotations thoroughly rinse gels
- Compatible with BCM1645 bungee cords to secure trays to platform
- Safe for incubator or cold room use



#### **Ordering Information**

Product Code	Description
BCM1618	Mini Blot 3D Rocker, fixed speed, small 10.5" x 7.5" platform with flat mat, 115 V

#### **Technical Data:**

<b>Speed</b>	12rpm	<b>Dimensions (WxDxH)</b>	10.5 x 7.5 x 8 inches 27 x 19 x 20 centimeters
<b>Tilt Angle</b>	±8°	<b>Platform Size</b>	10.5 x 7.5 inches 27 x 19 cm
<b>Load Capacity</b>	4lbs 1.8kg	<b>Operating Temperature</b>	+4°C to +65°C
<b>Weight</b>	4.4lbs/2kg	<b>Warranty</b>	2 Years
<b>Electrical</b>	115V, 50-60Hz, 50W		