



Product Information

Microplate Shaker

Product information for BCM1620:

Introduction

Microplate shaker. Speed can be adjusted. 3mm orbit for thorough mixing of microplate contents. 2 year warranty.

Features:

- Includes platform for up to 4 microplates
- Adjustable speed
- Safe for incubator or cold room use



Ordering Information

Product Code	Description
BCM1620	Shaker MP with 4 position microplate platform, 100-240V (US plug)

Technical Data:

Speed Range	200-1500rpm;	Dimensions (WxDxH)	11 x 10.7 x 4.5 inches 28 x 27 x 12 centimeters
Orbit	3mm	Timer	1min to 99hr/continuous
Load Capacity	4 microplates	Operating Temperature	+4°C to +65°C
Weight	12lbs/5.5kg	Warranty	2 Years
Electrical	100V, 50-60Hz, 250W		