



Product Information Orbital Blotter

Product information for BCM1621:

Introduction

Orbital rocker with flat mat. Speed can be adjusted.
2 year warranty.

Features:

- Adjustable speed
- Large work surface, non-slip rubber mat
- Safe for incubator or cold room use



Ordering Information

Product Code	Description
BCM1621	Blotter with non-slip rubber mat platform, 14" x 12", 115V

Technical Data:

Speed Range	3-70rpm	Dimensions (WxDxH)	14 x 12 x 8 inches 35 x 30 x 20 centimeters
Orbit	19mm	Platform Size	14 x 12 inches 35 x 30 cm
Load Capacity	3 x 2L, or 9lbs	Operating Temperature	+4°C to +65°C
Weight	7.5lbs/3.4kg	Warranty	2 Years
Electrical	115V, 50-60Hz, 50W		