



Product Information

Incubator Shakers

Product information for BCM1626/BCM1627/BCM1628:

Introduction

Shaking incubators with non-slip rubber mats, perfect for flasks, tube racks, and a wide variety of other containers. Available in multiple sizes to suit your needs. A larger shaking incubator with refrigeration capabilities is also offered. 2 year warranty.

Features:

- Hinged lids
- Seamless stainless steel chamber
- Digital temperature and timer controls

Ordering Information



Product Code	Description
BCM1626	Shaker Mini with non-slip rubber mat, 115V
BCM1627	Shaker 10L with non-slip rubber mat, 115V
BCM1628	Shaker 10LR with non-slip rubber mat, refrigerated shaking incubator, 115V

Technical Data:

	BCM1626	BCM1627	BCM1628
Temperature Range	Ambient+5°C to +70°C, 0.1°C incr.	Ambient+5°C to +60°C, 0.1°C incr.	Ambient-15°C to +70°C, 0.1°C incr.
Temperature Accuracy	±0.1°C		
Temperature Uniformity	±0.25°C		
Speed Range	30-300rpm, 1rpm incr.		
Shaking Orbit	19mm		
Platform Dimensions (WxD)	9.5 x 11.5 inches 23.5 x 19 cm	17.8 x 17.8 inches 45 x 45 cm	17.8 x 17.8 inches 45 x 45 cm
Exterior Dimensions (WxDxH)	11 x 16.2 x 13 inches 28 x 41 x 33 cm	20 x 26 x 19 inches 51 x 56 x 47 inches	20 x 31 x 22 inches 51 x 79 x 56 inches
Timer	1min to 48hr / Continuous	1min to 48hr / Continuous	1min to 48hr / Continuous
Weight	25 lbs / 11.5 kg	120 lbs / 55kg	170 lbs / 78kg
Electrical	115V, 50-60 Hz, 300W	120V, 50-60 Hz, 450W	120V, 50-60 Hz, 950W
Warranty	2 years		