



## **Product Information**

### **Strip Microcentrifuge**

#### **Product information for BCM1709/BCM1710/BCM1711:**

#### **Introduction**

Low profile mini centrifuge, optimized for 12 position PCR strips. 1 year warranty.

#### **Features:**

- Quiet motor
- Holds 4 x 12 position PCR strips, or 48 0.2ml PCR tubes
- Starts and stops with opening and closing of the lid
- Optional round rotor available for use with 1.5 and 2.0ml tubes



#### **Ordering Information**

Product Code	Description
BCM1709	12 Mini Centrifuge, holds 12 position PCR strips, 100-240V (US Plug)
BCM1710	Round rotor, 8 position
BCM1711	Adapters for 0.5ml tubes, pack of 8, for use with round rotor

#### **Technical Data:**

<b>Speed</b>	5,500rpm/2,000xg	<b>Dimensions (WxDxH)</b>	5.5 x 7.9 x 4.4 inches 14 x 20 x 11.2 centimeters
<b>Capacity</b>	48 x 0.2ml PCR Tubes 4 x PCR Strips (12 x 0.2ml)	<b>Weight</b>	2.2lbs/1kg
<b>Electrical</b>	Universal Voltage 100-240V, 50-60Hz, 300W	<b>Warranty</b>	1 Year