



# Product Information MicroPlate Microcentrifuge

## Product information for BCM1749/BCM1750:

### Introduction

Microplate mini centrifuge. Accepts all common formats of microplate (incl. skirted, semi and non). Concentrates reactants in bottom of wells for higher concentrations and optimized yields. Swing-out rotor eliminates the need to seal plates. 1 year warranty.

### Features:

- Swing-out rotor: Plates are angled while centrifuge is at rest but vertical while spinning
- Small footprint
- Eliminates the need for sealing film before centrifugation
- Decelerates in 7 seconds
- Optional tube adapter lets you spin PCR strips and tubes with the microplate microcentrifuge



Before & after centrifugation (angled position)



During centrifugation (horizontal position)

### Ordering Information

Product Code	Description
BCM1749	MicroPlate Microcentrifuge, with rotor and plate carriers, 115V
BCM1750	Tube adapter, holds 96 x 0.2ml or 8 and 12 position PCR strips

### Technical Data:

<b>Speed</b>	2,100rpm/400xg	<b>Dimensions (WxDxH)</b>	9.2 x 10.2 x 7.75 inches 23 x 26 x 19.7 centimeters
<b>Capacity</b>	2 x 2 PCR plates 2 x microtiter plates 24 x 0.2ml PCR strips	<b>Weight</b>	9lbs/4.1kg
<b>Electrical</b>	Universal Voltage 115-230V, 50-60Hz, 250W	<b>Warranty</b>	1 Year