



Product Information **Microtube and Plate Cooler**


Product information for BCM1834:

Introduction

Keeps temperature-sensitive samples cool on the bench. Place in freezer overnight before use. 2 year warranty.

Features:

- Maintains 0°C for 4 hours
- Ideal for 1.5 and 2.0 ml tubes, 0.2ml PCR tubes and strips, and 96 well PCR plates
- Ideal for PCR workstations



Ordering Information

Product Code	Description
BCM1834	Microtube and PCR Plate Cooler

Technical Data:

Temperature	0°C for up to 4 hours	Dimensions (WxDxH)	5.6 x 5.6 x 2.5 inches 14.3 x 14.3 x 6.4 centimeters
Capacity (top)	36 x 1.5-2ml Tubes	Weight	2.2lbs/1kg
Capacity (bottom)	96 x 0.2ml PCR Tubes 12 x PCR Strips (8 x 0.2ml) 1 x 96 well PCR plate 10 x 1.5-2.0ml tubes	Warranty	2 Year