



## Product Information

**Tip Volume: 100-1000uL**

### ***Product information for BT222/BT222-BS/BT222-NS/BT222-RS:***

#### **Description**

Biobasic Inc. made pipette tips by high purity virgin polypropylene and are designed for sensitive applications. They are compatible with robotic and manual pipetting systems and are tested for compatibility with BRAND, Eppendorf®, Gilson®/Rainin®, and Thermo Labsystems/Finnpipette® pipettes.

#### **Features:**

- Ideal for Precise, Accurate Pipetting: Tips are straight and feature clean, Uniform bores and hydrophobic surfaces to reduce sample retention and increase accuracy and reproducibility
- Available in variety of sizes.
- Autoclavable.

Item No.	Description	Package
BT222	Blue/bag	1000/bag
BT222-BS	Blue/bag/Sterile	1000/bag
BT222-NS	Natural/ bag/Sterile	1000/bag
BT222-RS	Natural/Rack/Sterile Certificated RNase, DNase-Free.	100/Rack



#### **Technical Specifications**

Item No.	BT222	BT222-BS	BT222-NS	BT222-RS
Package	1000/Bag	1000/Bag	1000/Bag	100/RACK
Volume	100-1000µl	100-1000µl	100-1000µl	100-1000µl
Color	Blue	Blue	Natural	Natural
Plastic materials	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Sterile	NO	YES	YES	YES
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Length	71 mm	71 mm	71 mm	71 mm
Inner diameter	7.6 mm	7.6 mm	7.6 mm	7.6 mm
Outside diameter	9.6 mm	9.6 mm	9.6 mm	9.6 mm

#### **Storage**

Store at Room temperature.