



## Product Information

### **BT620-N 1.5 ML GRADUATED CENTRIFUGE STERILE TUBE**

#### **Description:**

1.5ml Polypropylene Graduated Tubes with flat top caps.

Unlike polystyrene tube, Polypropylene tubes are chemical inert withstand centrifugation speed of 12,000g. Autoclavable. Tubes with flat top cap are easy to mark.

#### **Features:**

- Polypropylene tube
- 1.5ml
- With flat top cap
- Suitable for cell culture
- Autoclavable
- Storage RT

#### **Technical Specifications**

	BT620-N
Max. volume	1.5ml
Bottom shape	conical
Color	Natural
Max RCF	12,000 g
Plastic materials	Polypropylene
DNase-free	No
RNase-free	No
Withstand Temperatures	-80°C~121°C
Sterile	No
Dimensions	Outer Diameter: 10.7mm
	Length: 4.1mm
Qty./Pk	100X10