

Product Information

Graduated Plastic Cylinders

Description:

Graduated Plastic Cylinder is made with Polypropylene material. Polypropylene cylinders are chemical inert to most acids, bases and common solvents. They are autoclavable.



Features:

- Accuracy: The minimum division value is 1% of the full scale.
- The graduations: The graduations are slightly protuberant without paint.
- The foot: Non-detachable.
- Configuration of the foot: Hexagonal.
- Storage RT

Technical Specifications:

Cat.No.	CD001-1	CD001-2	CD001-3	CD001-4	CD001-5
Description	Graduated Plastic Cylinders	Graduated Plastic Cylinders	Graduated Plastic Cylinders	Graduated Plastic Cylinders	Graduated Plastic Cylinders
Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Capacity(ml)	100 ml	250 ml	500 ml	1 L	2 L
Color	Half Clear	Half Clear	Half Clear	Half Clear	Half Clear
Graduations	Yes	Yes	Yes	Yes	Yes
Cylinder Scale	10~100 ml	25~250 ml	50~500 ml	100~1000 ml	200~2000ml
Graduation Interval	1 ml	2.5 ml	5 ml	10 ml	20 ml
Tolerance +/-:	0.5 ml	1 ml	2 ml	5 ml	5 ml
Top Style	Funnel Top	Funnel Top	Funnel Top	Funnel Top	Funnel Top
Height	250 mm	340 mm	360 mm	440 mm	545 mm
Outer Diameter	32 mm	43 mm	54 mm	70 mm	90 mm