



Product Information

EZ-LINE Serological Pipettes, Volume: 1-25mL

Product information for CR4111/CR4112/CR4113/CR4114/CR4115:

Description

Sterile, individually-wrapped stripettes for tissue culture use. Made of clear polystyrene for easy visualization of fluid volume. Non-pyrogenic. Non-cytotoxic.

Features:

- Fine, discrete gradations
- Sterile; individually wrapped in paper/plastic wrap
- One-piece polystyrene construction.



Item No.	Description	Package
CR4111	EZ-LINE Serological Pipette, 1ml stripette, sterile, indiv. wrapped, 500 per case	500/case
CR4112	EZ-LINE Serological Pipette, 2ml stripette, sterile, indiv. wrapped, 500 per case	500/case
CR4113	EZ-LINE Serological Pipettes, 5ml stripette, Sterile, indiv. wrapped, 250 per case	250/case
CR4114	EZ-LINE Serological Pipettes, 10ml stripette, sterile, indiv. wrapped, 200 per case	200/case
CR4115	EZ-LINE Serological Pipettes, 25ml stripette, sterile, indiv. wrapped, 100 per case	100/case

Technical Specifications

Item No.	CR4111	CR4112	CR4113	CR4114	CR4115
Package	500/case	500/case	250/case	200/case	100/case
Volume	1mL	2mL	5mL	10mL	25mL
Color	Clear	Clear	Clear	Clear	Clear
Plastic materials	Polystyrene	Polystyrene	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes	Yes	Yes
Length	271mm	271mm	294mm	294mm	297mm

Storage

Store at Room temperature.

QF 24 Rev 0