



Bio Basic Inc.

A world leader in serving science

Data sheet

Product Name: Drying Reagent Silica Gel
Product Code: DSG110

Items		Results
Appearance		Blue or light blue vitreous particle
Particle Size		2-5mm
Absorbance	RH=20% %≥	---
	RH=35% %≥	---
	RH=50% %≥	20.0
Heating Loss %≤		5.0
Abrasion Ratio %≤		10.0
Percent of Pass (Particle Size) %≥		90
Color change	RH=20% %≥	---
	RH=35% %≥	---
	RH=50% %≥	Pinky
	Drying still restore primary colors	