

Product Information

Catalogue Number: MSD17
Product Name: pUC57
Size: 50µg
Concentration: 0.5µg/µl
Storage Buffer: 10mM Tris-HCl (pH 8.0), 1mM EDTA.
Storage: -20°C

Description:

The pUC57 plasmid is one of the most commonly used cloning vectors for molecular biology. The full-length pUC57 DNA is 2,710 bp in length. For more information, see the pUC57 plasmid DNA map below.

Features:

1. High purity plasmid, advanced purification technology.
2. More than 90 percent of the total DNA is supercoiled for easy transformations.
3. Many common restriction enzyme sites for easy cloning.

