



Product Information

96-Well 0.2ml PCR Plate

Product information for PCR-960-CH-S:

Description

A rigid top plate minimizes plate distortion, assures a dependable fit with the thermal cycler, and allows for leak-proof seal with Micro-Mats or Cap Strips. Each well has a capacity of 200ul. Sterilized individually. DNase and RNase-free.

Item No.	Color	Package
PCR-960-CH-S	Natural, Chimney-top, Sterile	10/pk

Technical Specifications

	96-Well 0.2ml PCR Plate
Max. volume	200µl
Plate Feature	Chimney
Well Bottom shape	Conical
Well Shape	Round
Color	Natural
Plastic materials	Polypropylene
Lids Included	No
Withstand Temperatures	-80°C~121°C
DNase-free	Yes
RNase-free	Yes
Sterile	Yes
Dimensions (W x D x H)	118.0 x 78.0 x 22.0 mm

Storage

Store at Room temperature.