



Bio Basic Inc.

A world leader in serving science

## **PRODUCT INFORMATION**

### **PSB Solution (PD8117)**

**Product name:** PBS Solution (Phosphate buffered saline),

**Catalog Number:** PD8117

#### **Product Description**

This PBS formulation is recommended for use in immunohistology protocols.

Each package contains 4 L of 10 concentrated PBS or 40 liters of 1 x solution. Each 1 x PBS solution contains 137 mM Sodium Chloride, 2.7mM Potassium Chloride and 10mM Phosphate Buffer, pH 7.4-7.6 at 25°C.

**Storage Temperature:** Room temperature.

#### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.