



## Product Information Parafilm

### **Product information for PF002:**

A versatile, moisture-resistant, thermoplastic flexible plastic for a wide variety of clinical and industrial applications, *Parafilm M*® film has been used with confidence for years in scientific research and medical settings.

It is a proven fast, economical covering or stopper for test tubes, beakers, vials, petri dishes, flasks and various instruments and objects with irregular surfaces.



Stock No	Roll Size	Packing per cs	Cs Wt.	Cubic Ft.
PF002	4" X 125"/roll	12 Rolls	18 lb	.91
PF003	2" X 250"/roll	24 Rolls	36 lb	.91