


Product Information

Floating Microtube Racks

Product information for R021:



Technical Specifications

Cat.No.	R021-N	R021-R	R021-Y	R021-G
Color	Natural	Red	Yellow	Green
Capacity	24 2.0 ml tubes & 8 0.5 ml tubes	24 2.0 ml tubes & 8 0.5 ml tubes	24 2.0 ml tubes & 8 0.5 ml tubes	24 2.0 ml tubes & 8 0.5 ml tubes
Applies to	0.5 or 2.0 ml Microcentrifuge Tubes	0.5 or 2.0 ml Microcentrifuge Tubes	0.5 or 2.0 ml Microcentrifuge Tubes	0.5 or 2.0 ml Microcentrifuge Tubes
Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Racks or Boxes Style	Tubes Floating	Tubes Floating	Tubes Floating	Tubes Floating
Format	24+8 Well	24+8 Well	24+8 Well	24+8 Well
Lids Included	No	No	No	No

Storage

Room temperature.