



## Product Information Aluminum Sealing Foil

### ***Product information for SMT4406:***

#### **Description**

This type of material is ideal for manual sealing during PCR work and also for high throughput applications. Its adhesive backing makes it easy to apply. Will resist temperatures from -86 °C to +105 °C. It is pierceable with a pipettor tip for easy access to sample. DMSO resistant.

It is recommended to use the Plate Roller (TM113) for ensuring a perfect bond, eliminating the danger of evaporation. No heat sealer needed.



| Code    | Description  | Qty/Pk |
|---------|--------------|--------|
| SMT4406 | Peeling foil | 100    |

#### **Storage**

Room temperature.