


Product Information

Serological Pipettes, Volume: 5-50 ml

Product information for SP102116/SP102121/SP102136/SP102161:

Description

Sterile, individually-wrapped stripettes for tissue culture use. Ascending & descending scales facilitate reading of both dispensing and remaining volume. Made of clear polystyrene for easy visualization of fluid volume. Non-pyrogenic. Non-cytotoxic.



Features:

- Exclusive color code, same colors as the international codes
- Filter plug inserted
- Quick & easy identification
- Sterilized
- Individual packaging

Item No.	Description	Package
SP102116	Serological Pipettes, 5ml stripette	400/case
SP102121	Serological Pipettes, 10ml stripette	400/case
SP102136	Serological Pipettes, 25ml stripette	200/case
SP102161	Serological Pipettes, 50ml stripette	160/case

Technical Specifications

Item No.	SP102116	SP102121	SP102136	SP102161
Units per case	400/case	400/case	200/case	160/case
Units per pack	100/pack	100/pack	50/pack	40/pack
Volume	5mL	10mL	25mL	50mL
Line Color	Blue	Orange	Red	Purple
Plastic materials	Polystyrene	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes	Yes
Length	338mm	338mm	343mm	343mm

Storage

Store at Room temperature.

QF 24 Rev 0