



Product Information

EZ-LINE Multiwell Tissue Culture Plates, Treated

Product information for SP41117, SP41123, SP41135, SP41207:

Description

Multiwell plates for tissue culture work. All plates are produced with optically clear, high quality crystalline polymers. Plates are treated for optimal cell attachment and growth and designed for reduced cross-contamination.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- TC treated
- High strength



Item No.	Description	Package
SP41117	EZ-LINE 6-well TC Plate, TC Treated sterile, indiv. wrapped, 50 per case	50/case
SP41123	EZ-LINE 12-well TC Plate, TC Treated sterile, indiv. wrapped, 50 per case	50/case
SP41135	EZ-LINE 24-well TC Plate, TC Treated Sterile, indiv. wrapped, 50 per case	50/case
SP41207	EZ-LINE 96-well TC Plate, TC Treated sterile, indiv. wrapped, 50 per case	50/case

Technical Specifications

Item No.	SP41117	SP41123	SP41135	SP41207
Package	50/case	50/case	50/case	50/case
Volume	3mL	2mL	1mL	200uL
Growth Area	9.60cm ²	3.80cm ²	1.90cm ²	0.32cm ²
Color	Clear	Clear	Clear	Clear
Plastic materials	Polystyrene	Polystyrene	Polystyrene	Polystyrene
Sterile	Yes	Yes	Yes	Yes
TC Treated	Yes	Yes	Yes	Yes

Storage

Room temperature.

QF 24 Rev 0