



## **Product Information** **EZ-LINE Cell Scraper**

### ***Product information for SP91131:***

#### **Description**

Bio Basic's Cell Scrapers are specially designed to facilitate collection of attached cells from a culture dish. Individual packing also eliminates possible chances of contamination during cell culture.

#### **Features:**

Sterilization by irradiation and individual packing eliminate possible chances of contamination during cell culture and scraper handling.



Item No.	Description	Package
SP91131	Cell Scraper	100/case

### **Technical Specifications**

Item No.	<b>SP91131</b>
Quantity per case	100/case
Packaging	Individually wrapped
Material (Handle / Blade)	PP / PE
Length (mm)	237
Blade Width (mm)	13
Sterile	Yes

#### **Storage**

Room temperature.



### Certificate of Compliance

**Product Description:** Cell culture products  
**Expiration Date:** 5 years from manufacture date

**Visual Attribute:** PASS  
**Functional Test:** PASS  
**Sterilization:** Sterilization was performed by exposed gamma irradiation according to ISO 11137 based on the *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization <SAL 10-6>*  
**Cytotoxicity Test:** Cytotoxicity test is conducted to qualify all material resins using USP class VI<87> and/or ISO 10993 standards for cytotoxicity and have been shown to be non-toxic  
**Pyrogens Test:** Non-pyrogenic tested according to the principles of the LAL test described in USP<85>/ 0.5EU/ml  
**RNase/DNase Test:** This product has been tested and is free of any detectable RNase/DNase contamination

- ✓ Bio Basic Inc. hereby certifies that the products identified above have been manufactured & inspected in compliance with products quality specifications and requirements as documented in our ISO 9001:2008 Quality Management System.