



Bio Basic Inc.

A world leader in serving science

Product Information

EZ-LINE Cell Scraper

Product information for SP91141:

Description

Designed for maximum contact with growth areas to easily scrape away attached cells. Rotatable blades for ease of use in round culture dishes.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- Rotatable blades
- Individually wrapped to minimized cross-contamination



Item No.	Description	Package
SP91141	EZ-LINE Cell Scraper, 2x30cm Sterile, indiv. wrapped, 100 per case	100/case

Technical Specifications

Item No.	SP91141
Quantity per case	100/case
Packaging	Individually Wrapped
Blade Width	2cm
Length	30cm
Color	Green handle, Clear blade
Plastic materials	Polypropylene handle, Polyethylene blade
Sterile	Yes

Storage

Room temperature.

QF 24 Rev 0

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com