




Product Information

EZ-LINE EZ-LINE Cell Lifter

Product information for SP91151:

Description

This cell lifter is designed for manual harvesting of attached cells from multiwell plates and culture dishes. The lifter has two blades of different sizes for use with different growth surfaces. The wide blade (23mm) is suitable for culture dishes and 6-well plates. The narrow blade is suitable for 12- and 24-well plates.



Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- Two blades for use with various growth surfaces
- Individually wrapped to minimized cross-contamination

Item No.	Description	Package
SP91151	EZ-LINE Cell Lifter, Double Blade, 23mm & 12mm Blades, 21.8cm Handle Sterile, indiv. wrapped, 100 per case	100/case

Technical Specifications

Item No.	SP91151
Quantity per case	100/case
Packaging	Individually Wrapped
Type	Double Blade
Blade Width	23mm and 12mm
Length	21.8cm
Color	Translucent White
Plastic materials	Polypropylene
Sterile	Yes

Storage

Room temperature.