



Product Information

TCF70-250 250mL Tissue Culture Bottle

Description:

Disposable. Optically clear polystyrene with flat distortion-free surface. Flasks have been tissue culture treated for optimum cell attachment. Graduated with screw cap. Sterile.

Features:

- Polystyrene
- 250mL
- Sterile
- Storage RT

Cat NO	Volume	Surface Area	Sterile or Not	Pack size
TCF70-250	250ml	75cm ²	Sterile	5/bag

Sterility testing:

Testing description	Accept value	Results
Sterility	pass	Pass
Free of endotoxins	Less than 1.1 x 10 ⁻¹² g/item	Pass
ATP-Free	Less than 1 x 10 ⁻¹⁵ g/item	Pass