



## Product Information

### TCF70-50 50mL Tissue Culture Bottle

#### **Description:**

Disposable. Optically clear polystyrene with flat distortion-free surface. Flasks have been tissue culture treated for optimum cell attachment. Graduated with screw cap. Sterile.

#### **Features:**

- Polystyrene
- 50mL
- Sterile
- Storage RT

Cat NO	Volume	Surface Area	Sterile or Not	Pack size
TCF70-50	50ml	25cm <sup>2</sup>	Sterile	10/bag

#### **Sterility testing:**

Testing description	Accept value	Results
Sterility	pass	Pass
Free of endotoxins	Less than 1.1 x 10 <sup>-12</sup> g/item	Pass
ATP-Free	Less than 1 x 10 <sup>-15</sup> g/item	Pass