



Product Information
5ml Polypropylene Round-bottomed Culture Tube

Product information for TD535 & TD535-S:

Description:

5ml Polypropylene Round-bottomed Culture Tube. Unlike polystyrene tube, Polypropylene tubes are chemical inert. Autoclavable. Tubes with screw cap are suitable for cell culture.

Features:

- Polypropylene tube
-5ml
-Dual cap
-Suitable for cell culture
-Autoclavable



Table with 4 columns: Item No., Description, Sterile or Not, Package. Rows include TD535 and TD535-S.

Storage

Room temperature.

Technical Specifications

Table with 3 columns: Specification, TD535, TD535-S. Rows include Tube Material, Tube Color, Capacity, Bottom Style, Self-Standing, Graduations, Graduation Range, Graduation Interval, Cap Material, Cap Style, Cap Color, Max RCF, Withstand Temperatures, Sterile, Dnase-free, Rnase-free, Outer Diameter, Length, Cap Diameter.