



Product Information

2ML SKIRTED TUBE WITH SCREW CAP, STERILE

Product information for TK361-TK367:

Description:

2.0ml Polypropylene Microtubes with screw caps. Unlike polystyrene tube, Polypropylene tubes are chemical inert withstand centrifugation speed of 12000xG. Autoclavable. Tubes with screw cap are suitable for cell culture.

Features:

- Polypropylene tube
- 2.0ml
- screw cap with O-ring
- Less than 12000xG for centrifugation speed
- Suitable for cell culture
- Autoclavable



Item No.	Description	Color/Sterile	Package
TK361-N	Skirted+ deep cap	natural, non sterile	500/bag
TK362-SN	Skirted+ deep cap	natural, sterile	500/bag
TK363-SB	Skirted+ deep cap	blue, sterile	500/bag
TK364-SR	Skirted+ deep cap	red, sterile	500/bag
TK365-SW	Skirted+ deep cap	white, sterile	500/bag
TK366-SY	Skirted+ deep cap	yellow, sterile	500/bag
TK367-SG	Skirted+ deep cap	green, sterile	500/bag

Storage

Room temperature.



Bio Basic Inc.

A world leader in serving science

Technical Specifications

	TK361-N	TK362-SN	TK363-SB	TK364-SR	TK365-SW	TK366-SY	TK367-SG
Tube Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Tube Color	Natural	Natural	Natural	Natural	Natural	Natural	Natural
Capacity(ml)	2ml	2ml	2ml	2ml	2ml	2ml	2ml
Bottom Style	Conical	Conical	Conical	Conical	Conical	Conical	Conical
Self-Standing	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Graduations	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Graduation Range	0.5~1.5ml	0.5~1.5ml	0.5~1.5ml	0.5~1.5ml	0.5~1.5ml	0.5~1.5ml	0.5~1.5ml
Graduation Interval	0.25ml	0.25ml	0.25ml	0.25ml	0.25ml	0.25ml	0.25ml
Cap Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Cap Style	screw cap with O-ring	screw cap with O-ring	screw cap with O-ring	screw cap with O-ring	screw cap with O-ring	screw cap with O-ring	screw cap with O-ring
Cap Color	Natural	Natural	Blue	Red	White	Yellow	Green
Max RCF	12000×G	12000×G	12000×G	12000×G	12000×G	12000×G	12000×G
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Sterile	No	Yes	Yes	Yes	Yes	Yes	Yes
Dnase-free	No	Yes	Yes	Yes	Yes	Yes	Yes
Rnase-free	No	Yes	Yes	Yes	Yes	Yes	Yes
Silicone Washer	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Cap Diameter	12mm	12mm	12mm	12mm	12mm	12mm	12mm

OF 24 Rev 0

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
 Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
 Email: order@biobasic.com Web: www.biobasic.com