



PRODUCT INFORMATION

5 ml Polypropylene Round-bottom Culture Tubes

Product information for TP536:

Description:

5ml Polypropylene Round-bottom Culture Tubes with Snap Caps. Unlike polystyrene tube, Polypropylene tubes are chemical inert. Autoclavable. Suitable for cell culture.

Item No.	Description	Sterile or Not	Package
TP536-N	Tube & natural cap	Non sterile	100/bag
TP536-SN	Tube & natural cap	Sterile	50/bag
TP536-B	Tube & blue cap	Non sterile	100/bag
TP536-SB	Tube & natural cap	Sterile	50/bag
TP536-R	Tube & red cap	Non sterile	100/bag
TP536-SR	Tube & natural cap	Sterile	50/bag
TP536-Y	Tube & yellow cap	Non sterile	100/bag
TP536-SY	Tube & natural cap	Sterile	50/bag



Storage: Room temperature.

Technical Specifications

	TP536-N	TP536-SN	TP536-B	TP536-SB	TP536-R	TP536-SR	TP536-Y	TP536-SY
Tube Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Tube Color	Natural	Natural	Natural	Natural	Natural	Natural	Natural	Natural
Capacity	5ml	5ml	5ml	5ml	5ml	5ml	5ml	5ml
Graduations	No	No	No	No	No	No	No	No
Tube Style	Culture	Culture	Culture	Culture	Culture	Culture	Culture	Culture
Bottom Style	Round	Round	Round	Round	Round	Round	Round	Round
Self-Standing	No	No	No	No	No	No	No	No
Neck Style	Cylinder	Cylinder	Cylinder	Cylinder	Cylinder	Cylinder	Cylinder	Cylinder
Top Style	Snap Cap	Snap Cap	Snap Cap	Snap Cap	Snap Cap	Snap Cap	Snap Cap	Snap Cap
Cap Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Cap Style	Plug Seal	Plug Seal	Plug Seal	Plug Seal	Plug Seal	Plug Seal	Plug Seal	Plug Seal
Cap Color	Natural	Natural	Blue	Blue	Red	Red	Yellow	Yellow
Max RCF	-	-	-	-	-	-	-	-
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Sterile	No	Yes	No	Yes	No	Yes	No	Yes
Dnase-free	No	Yes	No	Yes	No	Yes	No	Yes
Rnase-free	No	Yes	No	Yes	No	Yes	No	Yes
Outer Diameter	11.9mm	11.9mm	11.9mm	11.9mm	11.9mm	11.9mm	11.9mm	11.9mm
Length	85mm	85mm	85mm	85mm	85mm	85mm	85mm	85mm
Cap Diameter	12.10mm	12.10mm	12.10mm	12.10mm	12.10mm	12.10mm	12.10mm	12.10mm