



CERTIFICATE OF ANALYSIS

Product	Acrylamide/Bisacrylamide (37.5:1) Solution 40% (w/v)
Grade	Ultra Pure Grade
Product Code	A0008
Lot	

Test Items	Specifications	Actual Results
Appearance	Clear, colorless solution	
Concentration of acrylamide (w/v)	38.96 %	
Concentration of bisacrylamide (w/v)	1.04 %	
Conductivity (@18°C)	≤ 5 μmho	

Storage: 4 °C, darkness.