



CERTIFICATE OF ANALYSIS

Product	TBS (TRIS-BUFFER SALINE) 10X PREMIXED POWDER
Grade	BIOTECH GRADE
Product Code	A0025
Formula	N/A
MW	N/A
Lot No	

Test Items	Specifications	Results
Appearance	White powder	
Conductivity(1.2%,Water)	REPORT	
pH @25C(1.2%, Water)	7.25-7.55	
Solubility (1.2%, Water)	PASS	

Note: Each pack makes 2L of 10X concentrate, or 20L of 1X TBS buffer.

Storage: Store at room temperature. Protect from moisture.