



---

## CERTIFICATE OF ANALYSIS

---

Product	TG-SDS Buffer (TRIS-GLYCINE-SDS) 10x Solution
Grade	N/A
Product Code	A0030
Formula	N/A
MW	N/A
Lot No	

---

Test Items	Specifications	Results
Appearance	Liquid	
pH (1 x solution)	8.3~8.7	
Usage	Single usage only	

---

Description: A single strength (1X) solution contains 0.025M Tris base, 0.192M Glycine and 0.1% Sodium Dodecyl Sulfate.

Storage: R.T.