



CERTIFICATE OF ANALYSIS

Product	TG-SDS Buffer (Trisbase-Glycine-SDS) 10X Premix powder
Grade	Biotech
Product Code	A0031
Formula	N/A
MW	N/A
Lot No	

Test Items	Specifications	Results
Appearance	White free-flowing powder	
pH (5X Dilution of Reconstituted Buffer)	8.1-8.6	

Description: A single strength (1X) solution contains 0.025M Tris base, 0.192M Glycine and 0.1% sodium dodecyl sulfate. Single usage only.

Note: 1 pack makes 4 liters of 10X concentrate TG-SDS buffer or 40 liters of 1X TG-SDS buffer

Storage: R.T.